

Calcimeter



Brand: OFI Testing Equipment, Inc.

Product Code: 152-95

Availability: 6

Lead Time: 120

Description

The Calcimeter accurately and quickly determines if scale build up is composed of calcium carbonate. The calcite to dolomite content of the unknown sample aids in determining which chemical treating program to implement. The entire test procedure requires 15 to 30 minutes for both calcite and dolomite determination.

Part number: #152-95

Size: 13" × 8" × 2.5" (33 × 20 x 6 cm)

Weight: 2 lb 7 oz (1.1 kg)