

# Visual Examination of Used Electrical Insulating Oils

## Visual Examination of Used Electrical Insulating Oils of Petroleum Origin in the Field

### test method

Provides an estimate of the color and condition of in-service oils by visual observation and comparison with ASTM color standards in an oil comparator.

### oil comparator

- Conforms to ASTM D1524 specifications
- Yields results equivalent to ASTM D1500

Complete ASTM oil color test outfit for comparison of oils against ASTM color standards. Includes two color discs, ranging from 0.5 to 5.0 in 10 steps and 5.0 to 8.0 in 7 steps. Magnifying prism brings the sample and standard color fields together for side by side comparison. Portable unit is suitable for laboratory or field use. Supplied with two precision 33mm rectangular glass cells, carrying case and instructions.

### specifications

#### Shipping Information

Shipping Weight: 10 lbs (4.5kg)

Dimensions: 1 Cu. ft.



K13203 Oil Comparator

### ordering information

catalog no.	description	qty
K13203	Oil Comparator	1
<b>accessories</b>		
K13204	Daylight Illuminator, 115V Provides uniform lighting for Oil Comparator	1
K13294	Daylight Illuminator, 220-240V	
K13205	Rectangular Glass Cell	