


PowderSafe™ Type A Enclosure

110V



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene - No Rust
	Walls	Polycarbonate - No Rust
	Front Sash	Clear Polycarbonate
	Head Assembly / Filtration Housing	Molded ABS
	Initial HEPA Filter	Installed
Controls		ON/OFF
Face Velocity - Variable		2-Speed, High-Low; Door Controlled
Lighting		Vapor-Proof Fluorescent
Electrical Cord Pass-Through		Two; Right & Left Side Ports
Blowers		Brushless-Sparkless
HEPA Filter Monitoring Package		Audible / Visual Alarm

DIMENSIONS												
Product #	External Dimensions					Internal Dimensions			Sash Opening Dimensions			
	Width at Base	Max Width	Depth at Base	Max Depth	Height	Width	Depth	Height	Width	Height		
										Opening	Sash Fold 1	Sash Fold 2
AC705	23 3/4"	24 1/2"	24"	25 1/4"	29 1/2"	23 1/4"	23 1/2"	19 1/2"	21 1/4"	4"	9 1/2"	19"

TECHNICAL SPECIFICATIONS							
Product #	Electrical			Volume of Filtered Air	Prefilters	Filtration Layer	Unit Weight*
	Voltage	Amps	Power				
AC705	110V	1.0 A	110 W	175 CFM	1	1 HEPA	110 lbs.

**Unit weight is an approximation*



AirClean[®]Systems
The Fume Control Experts.



www.AirCleanSystems.com • info@aircleansystems.com
Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015