

**Dead Air Box**

**110V**



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene - No Rust
	Walls	Clear Polycarbonate
	Front Sash	Clear Polycarbonate
	Pipette Shelf	Clear Polycarbonate
	Head Assembly	Molded ABS
Controls		UVTect™ Microprocessor Controller
Lighting		Built-in Fluorescent
UV Lighting		254-nm Shortwave UV
UV Monitoring Package		Audible / Visible Alarm
UV Safety Switch		Included
Automatic UV Timer		Custom Timed Irradiation
Lab Event Timer		Adjustable Audible / Visible Timer

DIMENSIONS											
Product #	External Dimensions				Internal Dimensions			Sash Opening Dimensions			
	Width at Base	Max Width	Depth at Base	Height	Width	Depth	Height	Width	Height		
									Sash Fold 1	Sash Fold 2	
AC632DBC	31 3/4"	32 1/2"	24"	29 1/2"	31 1/4"	23 1/2"	21 1/2"	29 1/2"	9 1/2"	19 1/2"	
AC632TDBC	31 3/4"	32 1/2"	24"	39 1/2"	31 1/4"	23 1/2"	30 1/2"	29 1/2"	9 1/2"	25 3/4"	
AC648DBC	47 3/4"	48 1/2"	24"	29 1/2"	47 1/4"	23 1/2"	21 1/2"	45 1/4"	9 1/2"	19 1/2"	
AC648TDBC	47 3/4"	48 1/2"	24"	39 1/2"	47 1/4"	23 1/2"	30 1/2"	45 1/4"	9 1/2"	25 3/4"	

TECHNICAL SPECIFICATIONS				
Product #	Electrical			Unit Weight*
	Voltage		Power	
AC632DBC	110V		80 W	100 lbs.
AC632TDBC	110V		80 W	115 lbs.
AC648DBC	110V		140 W	130 lbs.
AC648TDBC	110V		140 W	145 lbs.

\*Unit weight is an approximation



**www.AirCleanSystems.com • info@aircleansystems.com**  
**Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015**