

MailSafe™ Ductless Class 1 Mail Handling Workstation 110V



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene - No Rust
	Walls	Clear Polycarbonate
	Front Sash	Clear Polycarbonate
	Head Assembly	Molded ABS
HEPA Filters		Initially Included; 99.99% Efficient at 0.3 µm
Controls		AirSafe™ Automatic Safety Controller
Lighting		Built-in Fluorescent
Blowers		Brushless-Sparkless; Automatically Engaged
Lab Event Timer		Adjustable Audible / Visual Timer
HEPA Filter Monitoring Package		Visual Helic Gauge Indicator
HEPA Filter Change Alarm		Audible / Visual Alarm

DIMENSIONS											
Product #	External Dimensions				Internal Dimensions			Sash Opening Dimensions			
	Width at Base	Max Width	Depth	Height	Width	Depth	Height	Width	Height		
									Sash Fold	Sash Fold	
AC632MS	31 ¾"	32 ½"	24"	29 ½"	31 ¼"	23 ½"	19 ½"	29 ¼"	9 ½"	19"	
AC648MS	47 ¾"	48 ½"	24"	29 ½"	47 ¼"	23 ½"	19 ½"	45 ¼"	9 ½"	19"	

TECHNICAL SPECIFICATIONS						
Product #	Electrical			Volume of Filtered Air	Filtration Layer	Unit Weight*
	Voltage	Amps	Power			
AC632MS	110V	2.0 A	220 W	480 CFM	1 HEPA	150 lbs.
AC648MS	110V	4.1 A	450 W	960 CFM	2 HEPAs	230 lbs.

*Unit weight is an approximation



www.AirCleanSystems.com • info@aircleansystems.com
Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015