

Polypropylene Total Exhaust Fume Hood

110V



| GENERAL SPECIFICATIONS | | |
|------------------------|-------------------|---|
| Construction Material | Work Surface | Polypropylene - No Rust |
| | Walls | Polypropylene - No Rust |
| | Front Sash | Clear Polycarbonate |
| | Metal Free Screws | TPU |
| Roof Blower | | Optional, Variable-Speed Polypropylene Blower |
| Lighting | | Vapor-Proof Fluorescent |

| DIMENSIONS | | | | | | | | | | |
|------------|---------------------|-------|--------|---------------------|-------|--------|-------------------------|---------|-------------|-------------|
| Product # | External Dimensions | | | Internal Dimensions | | | Sash Opening Dimensions | | | |
| | Width | Depth | Height | Width | Depth | Height | Width | Opening | Height | |
| | | | | | | | | | Sash Fold 1 | Sash Fold 2 |
| AC3000TE | 36" | 36" | 59" | 35" | 32" | 43" | 31" | 10" | 18" | 31" |
| AC4000TE | 48" | 36" | 59" | 47" | 32" | 43" | 43" | 10" | 18" | 31" |
| AC5000TE | 60" | 36" | 59" | 59" | 32" | 43" | 55" | 10" | 18" | 31" |
| AC6000TE | 72" | 36" | 59" | 71" | 32" | 43" | 67" | 10" | 18" | 31" |

| TECHNICAL SPECIFICATIONS | | | | | |
|--------------------------|------------|-------|-------|-----------------------------------|--------------|
| Product # | Electrical | | | Volume of Air Required at 100 LFM | Unit Weight* |
| | Voltage | Amps | Power | | |
| AC3000TE | 110V | 0.4 A | 50 W | 225 CFM | 350 lbs. |
| AC4000TE | 110V | 0.6 A | 60 W | 310 CFM | 450 lbs. |
| AC5000TE | 110V | 0.6 A | 70 W | 400 CFM | 570 lbs. |
| AC6000TE | 110V | 0.6 A | 70 W | 475 CFM | 700 lbs. |

* Unit weight is an approximation



AirClean® Systems



www.aircleansystems.com • info@aircleansystems.com

Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015