

Mini CyanoSafe™ Cyanoacrylate Fuming Chamber

110V



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene
	Chamber	Polypropylene
	Door	Polycarbonate with Handle Lock
	Hanging Rods & Clips	Stainless Steel
	Filtration Housing	Polypropylene
Controls		AirSafe™ Cyano Safety Controller
Fuming Cycle		Automatic Processing with Safety Switch Start
Lighting		Vapor-Proof Fluorescent
Blowers		Brushless-Sparkless
Humidity Sensor		Included
Hot Plate		Included
Carbon Filter Monitoring Package		Audible / Visual Alarm
Test Strip Holder		Included
Certifications		UL and CSA

DIMENSIONS							
Product #	External Dimensions				Internal Dimensions		
	Width	Depth at Base	Max Depth	Height	Width	Depth	Height
CAS24	24"	17 3/8"	18 1/2"	32"	23 1/2"	16 3/4"	17"

TECHNICAL SPECIFICATIONS						
Product #	Electrical			Filtration Layer	Unit Weight*	
	Voltage	Amps	Power			
CAS24	110V	6.2 A	680 W	1 Carbon	120 lbs.	

*Unit weight is an approximation



AirClean® Systems
The Fume Control Experts.



www.AirCleanSystems.com • info@aircleansystems.com
Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015