

AC600 Series Ductless Endoscopy Workstation

110V



GENERAL SPECIFICATIONS		
Construction Material	Spill Tray	Welded Polypropylene Spill Tray
	Walls	Clear Polycarbonate
	Front Sash	Clear Polycarbonate
	Head Assembly / Filtration Housing	Molded ABS
	Carbon Filter	Initially Included; Glutaraldehyde / OPA Filter
	Mobile Cart	Optional; Epoxy Powder-Coated Stainless Steel Cart
Controls		AirSafe™ Automatic Safety Controller
Face Velocity - Variable		User Defined; 40-120 LFM
Lighting		Vapor-Proof Fluorescent*
Blowers		Variable-Speed, Brushless-Sparkless
Low Airflow Alarm		Custom Audible / Visible Alarm
Peristaltic Pump		Optional
Soaking Cycle Timer		Adjustable Audible / Visual Timer
Filter Monitoring Package		Audible / Visible Alarm
Certifications		UL and CSA

DIMENSIONS										
Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions			
	Width	Depth	Height	Width	Depth	Height	Width	Height		
								Opening	Fold 1	Fold 2
AC632C	32"	24"	29 ½"	31 ¼"	23 ½"	19 ½"	29 ½"	4"	9 ½"	19"
AC648C	48"	24"	29 ½"	47 ¼"	23 ½"	19 ½"	45 ½"	4"	9 ½"	19"

TECHNICAL SPECIFICATIONS							
Product #	Electrical			Volume of Filtered Air	Prefilters Required	Primary Filtration Layer	Unit Weight**
	Voltage	Amps	Power				
AC632C	110V	1.0 A	100 W	175 CFM	1	1 Carbon	180 lbs.
AC648C	110V	2.0 A	200 W	350 CFM	2	2 Carbons	260 lbs.

*Fluorescent lighting is optional on 32" models and will add 2 ¾" external depth

**Unit weight is an approximation



www.AirCleanSystems.com • info@aircleansystems.com

Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015