

AC600 Series HEPA-Filtered Workstation with UV Light

110V



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene - No Rust
	Walls	Polycarbonate - No Rust
	Front Sash	Clear Polycarbonate
	Head Assembly	Molded ABS
Controls		UVTect™ Microprocessor Controller
Lighting		Built-In Fluorescent
Blowers		Brushless-Sparkless; Automatically Engaged
Lab Event Timer		Adjustable Audible / Visual Timer
UV Lighting		254 nm Shortwave UV
Automatic UV Timer		Custom Timed Irradiation
UV Safety Switch		Included
UV Monitoring Package		Audible Visual Alarm
Initial HEPA Filter(s)		Installed; 99.99% Efficient at 0.3µm
HEPA Filter Monitoring Package		Audible / Visible Alarm

DIMENSIONS											
Product #	External Dimensions				Internal Dimensions			Sash Opening Dimensions			
	Width at Base	Max Width	Depth at Base	Height	Width	Depth	Height	Width	Height		
									Sash Fold 1	Sash Fold 2	
AC632UPUV	31 ¾"	32 ½"	24"	29 ½"	31 ¼"	23 ½"	19 ½"	29 ¼"	9 ½"	19 ½"	
AC648UPUV	47 ¾"	48 ½"	24"	29 ½"	47 ¼"	23 ½"	19 ½"	45 ¼"	9 ½"	19 ½"	

TECHNICAL SPECIFICATIONS						
Product #	Electrical		Volume of Filtered Air	Prefilters	Filtration Layer	Unit Weight*
	Voltage	Power				
AC632UPUV	110V	135 W	175 CFM	1	1 HEPA	160 lbs.
AC648UPUV	110V	230 W	350 CFM	2	2 HEPA	260 lbs.

*Unit weight is an approximation



AirClean[®]Systems
The Fume Control Experts.



www.AirCleanSystems.com • info@aircleansystems.com
Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015