

CleanShield® Endoscope Storage Cabinet 6-Scope Storage Solution

110V



GENERAL SPECIFICATIONS		
Construction Material	Shell	Polypropylene - No Rust
	Door	Epoxy Powder-Coated Aluminum with Clear Polycarbonate Viewing Window
	Filtration Housing	Polypropylene - No Rust
Blowers		Brushless-Sparkless
Blower Speed Settings		Automatic High and Low
Controls		ON/OFF
Initial HEPA Filters		Installed
HEPA Filter Monitoring Package		Audible / Visual Alarm
Certifications		UL and CSA

DIMENSIONS							
Product #	External Dimensions			Internal Dimensions			
	Width	Depth	Height	Width	Depth	Height	Max. Shaft Length
ACCSES06	39"	15 1/2"	94"	30"	13"	89"	81 1/4"

TECHNICAL SPECIFICATIONS								
Product #	Electrical			Number of Scopes Stored	Prefilters	Filtration Layer	Unit Weight*	Shipping Weight*
	Voltage	Amps	Power					
ACCSES06	110 V	1.2 A	130 W	6	1	1 HEPA	150 lbs	290 lbs

*Weight is an approximation.



www.AirCleanSystems.com • info@aircleansystems.com
 Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015