

CleanShield® Endoscope Storage Cabinet 9-Scope Storage Solution

110V



GENERAL SPECIFICATIONS		
Construction Material	Shell	Polypropylene - No Rust
	Door	Epoxy Powder-Coated Aluminum with Clear Polycarbonate Viewing Window
	Filtration Housing	Polypropylene - No Rust
Blowers		Brushless-Sparkless
Blower Speed Settings		Automatic High and Low
Controls		ON/OFF
Initial HEPA Filters		Installed
HEPA Filter Monitoring Package		Audible / Visual Alarm
Certifications		UL and CSA

DIMENSIONS							
Product #	External Dimensions			Internal Dimensions			
	Width	Depth	Height	Width	Depth	Height	Max. Shaft Length
ACCSES09	39"	15 ½"	94"	10" (per section)	13"	89"	81 ¼"

TECHNICAL SPECIFICATIONS								
Product #	Electrical			Number of Scopes Stored	Prefilters	Filtration Layer	Unit Weight*	Shipping Weight*
	Voltage	Amps	Power					
ACCSES09	110 V	1.2 A	130 W	9	1	1 HEPA	175 lbs	315 lbs

*Unit weight is an approximation.



www.AirCleanSystems.com • info@aircleansystems.com
 Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015