

# Neutralizer Station

# 110V



GENERAL SPECIFICATIONS		
Construction Material	Lid	Polypropylene
	Stand	Epoxy Powder-Coated Steel
	Tank	Polypropylene
Controls		ON/OFF
Neutralizer Addition Indicator		Visual Indicator
Discharge Button		Visual Indicator
Initial Vapor-Management Filter		Included
Certifications		UL

DIMENSIONS			
Product #	External Dimensions		
	Width	Depth	Height
ACDS03	13 ½"	16"	41 ½"

TECHNICAL SPECIFICATIONS						
Product #	Electrical			Filtration Layer	Tank Volume	Unit Weight*
	Voltage	Amps	Power			
ACDS03	110V	0.5 A	50 W	1 Carbon	6.5 Gallons	80 lbs.

\*Unit weight is an approximation.



**AirClean<sup>®</sup>Systems**  
The Fume Control Experts.



[www.AirCleanSystems.com](http://www.AirCleanSystems.com) • [info@aircleansystems.com](mailto:info@aircleansystems.com)  
Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015