

CleanShield[®] TEE Probe Storage Cabinet

110V



GENERAL SPECIFICATIONS		
Construction Material	Shell	Polypropylene - No Rust
	Door	Polypropylene with Clear Polycarbonate Viewing Window
	Filtration Housing	Polypropylene - No Rust
	Lower Back Wall	Closed Cell Foam
Blowers		Brushless-Sparkless
Controls		ON/OFF
Initial HEPA Filters		Installed
HEPA Filter Monitoring Package		Audible / Visual Alarm
Certifications		UL and CSA

DIMENSIONS								
Product #	External Dimensions			Internal Dimensions				
	Width	Depth	Height	Width	Depth	Height	Height Shelf 1	Height Shelf 2
ACTE03	25 ¼"	10 ⅞"	88"	24 ½"	9"	70 ½"	41 ½"	53 ¾"
ACTE06	51 ½"	10 ⅞"	88"	49"	9"	70 ½"	41 ½"	53 ¾"

TECHNICAL SPECIFICATIONS							
Product #	Electrical		Number of Probes Stored	Prefilters	Filtration Layer	Unit Weight*	Shipping Weight*
	Voltage	Power					
ACTE03	110 V	100 W	3	1	1 HEPA	125 lbs	217 lbs
ACTE06	110 V	170 W	6	2	2 HEPA	250 lbs	217 lbs (per crate)

*Unit weight is an approximation.



www.AirCleanSystems.com • info@aircleansystems.com

Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015