

CleanShield® Endocavity Ultrasound Storage Cabinet

110V



GENERAL SPECIFICATIONS		
Construction Material	Shell	Polypropylene - No Rust
	Door	Polypropylene with Clear Polycarbonate Viewing Window
	Filtration Housing	Polypropylene - No Rust
Blowers		Brushless-Sparkless
Controls		ON/OFF
Initial HEPA Filters		Installed
HEPA Filter Monitoring Package		Audible / Visual Alarm
Certifications		UL and CSA

DIMENSIONS									
Product #	External Dimensions			Internal Dimensions					
	Width	Depth	Height	Width	Depth	Height	Height Shelf 1	Height Shelf 2	Height Shelf 3
ACVR06	25 3/4"	10 7/8"	88"	24 1/2"	9"	70 1/2"	39"	51 1/2"	53 1/4"

TECHNICAL SPECIFICATIONS								
Product #	Electrical			Number of Probes Stored	Prefilters	Filtration Layer	Unit Weight*	Shipping Weight*
	Voltage	Amps	Power					
ACVR06	110 V	0.9 A	100 W	6	1	1 HEPA	125 lbs	217 lbs

*Weight is an approximation.



www.AirCleanSystems.com • info@aircleansystems.com
 Phone: (919) 255-3220 • Toll Free: (800) 849-0472 • Fax: (919) 528-0015