

# Safety Data Sheet

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

## SECTION 1: Identification

### Product identifier

Trade name/designation:	Ammonium iron (III) citrate redbrown
Product No.:	1980
Synonyms:	none

### Relevant identified uses of the substance or mixture and uses advised against

<b>Recommended use</b>	For Laboratory, Research or Manufacturing Use.
<b>Uses advised against</b>	Not determined.

### Details of the supplier of the safety data sheet

#### Supplier

##### Avantor Performance Materials, LLC.

Street	100 Matsonford Rd, Suite 200
Postal code/City	Radnor, PA 19087, United States
Telephone	+1-855-282-6867
Telefax	+1-610-573-2610

#### Emergency phone number

Telephone	+1-800-424-9300 (Chemtrec, 24 hrs/day, 7 days/week, USA and Canada)
-----------	---

#### Preparation Information

Product Information Compliance

E-mail	SDS@avantorsciences.com
--------	-------------------------

## SECTION 2: Hazard identification

### Classification of the substance or mixture

This substance is classified as not hazardous according to regulation 29 CFR 1910.1200 (OSHA HCS).

### Label elements

According to regulation 29 CFR 1910.1200 (OSHA HCS) the product does not have to be labelled.

### Hazard(s) not otherwise classified (HNOC)

none

### SECTION 3: Composition/information on ingredients

#### Substances

Substance name	Ammonium iron (III) citrate redbrown
Molecular formula	$C_6H_{(5+4y)}Fe(x)N(y)O_7$
Molecular weight	no data available
CAS No.	1185-57-5

### SECTION 4: First aid measures

#### General information

When in doubt or if symptoms are observed, get medical advice. Change contaminated, saturated clothing. Wash contaminated clothing before reuse. Do not leave affected person unattended.

#### In case of inhalation

Remove casualty to fresh air and keep warm and at rest. Obtain medical attention if symptoms appear.

#### In case of skin contact

Gently wash with plenty of soap and water. In case of skin reactions, consult a physician.

#### After eye contact:

Rinse immediately carefully and thoroughly with eye-bath or water. Obtain medical attention if symptoms appear.

#### In case of ingestion

Rinse mouth thoroughly with water. Call a doctor if you feel unwell.

#### Most important symptoms/effects, acute and delayed

No known symptoms to date.

#### Indication of any immediate medical attention and special treatment needed

No special information on medical attention and special treatment available.

### SECTION 5: Fire fighting measures

#### Extinguishing media

##### Suitable extinguishing media

ABC-powder  
Carbon dioxide (CO<sub>2</sub>).  
Dry sand  
Nitrogen

#### Extinguishing media which must not be used for safety reasons

Water spray.  
Full water jet

#### Specific hazards arising from the chemical

Combustible substance.  
This material is combustible, but will not ignite readily.  
Fire may produce irritating, corrosive and/or toxic gases.  
In case of fire may be liberated:  
Pyrolysis products, toxic

#### Advice for firefighters

In case of fire and/or explosion do not breathe fumes.  
Protective equipment and precautions for firefighters:

Wear a self-contained breathing apparatus and chemical protective clothing.  
Co-ordinate fire-fighting measures to the fire surroundings.  
In case of fire: Evacuate area.

## SECTION 6: Accidental release measures

### Personal precautions, protective equipment and emergency procedures

For non-emergency personnel: Wear personal protection equipment (refer to section 8). Do not breathe gas/fume/vapor/spray. Avoid contact with skin, eyes and clothes. Stop leak if safe to do so. First Aid, decontamination, treatment of symptoms. For emergency responders: This material is combustible and can be ignited by heat, sparks, flames, or other sources of ignition (e.g. static electricity, pilot lights, or mechanical/electrical equipment). Wear fire/flame resistant/retardant clothing. Wear full chemical protective clothing. In case of major fire and large quantities: Evacuate area. Fight fire remotely due to the risk of explosion.

### Environmental precautions

No special measures are necessary.

### Methods and material for containment and cleaning up

Take up mechanically, placing in appropriate containers for disposal. Dispose according to legislation.

### Additional information

Personal protection equipment (PPE): see section 8 SECTION 13. Information regarding the disposal of the products

## SECTION 7: Handling and storage

### Precautions for safe handling

Advices on safe handling

No special measures are necessary.

Measures to prevent fire, aerosol and dust generation

No special measures are necessary.

Measures required to protect the environment

No special measures are necessary.

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

### Conditions for safe storage, including any incompatibilities

Recommended storage temperature: no data available

Storage: Store in a well-ventilated place. Store in a closed container. Keep only in original container.

Packaging materials: Polyethylene Unsuitable materials and coatings of containers/equipment: No information available.

## SECTION 8: Exposure controls/personal protection

### Control parameters

Does not contain substances above concentration limits fixing an occupational exposure limit.

## Engineering controls

### Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment. If handled uncovered, arrangements with local exhaust ventilation have to be used.

### Personal protection equipment (PPE)

Wear suitable protective clothing. When handling with chemical substances, protective clothing must be worn.

#### *Eye/face protection*

Eye glasses with side protection

#### *Skin protection*

Wear suitable gloves. When handling with chemical substances, protective gloves must be worn. In the case of wanting to use the gloves again, clean them before taking off and air them well. Check leak tightness/impermeability prior to use.

#### By short-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,12 mm
Breakthrough time	480 min

#### By long-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,12 mm
Breakthrough time	480 min

#### *Respiratory protection*

no data available

#### *Additional information*

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

#### *Environmental exposure controls*

no data available

## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

(a) Appearance	
Physical state:	solid
Color:	brown
(b) Odor:	characteristic
(c) Odor threshold:	no data available

### Safety relevant basic data

(d) pH:	6-8
(e) Melting point/freezing point:	no data available
(f) Initial boiling point and boiling range:	no data available
(g) Flash point:	no data available
(h) Evaporation rate:	no data available
(i) Flammability (solid, gas):	not applicable
(j) Flammability or explosive limits	
Lower explosion limit:	no data available
Upper explosion limit:	no data available
(k) Vapor pressure:	no data available
(l) Vapor density:	no data available
(m) Density:	0.73 g/cm <sup>3</sup> (20 °C)
(n) Solubility(ies)	
Water solubility:	1200 g/l (20°C)
Soluble (g/L) in Ethanol:	no data available
(o) Partition coefficient: n-octanol/water:	no data available
(p) Auto-ignition temperature:	510 °C
(q) Decomposition temperature:	not applicable
(r) Viscosity	
Kinematic viscosity:	no data available
Dynamic viscosity:	no data available
(s) Explosive properties:	not applicable
(t) Oxidising properties:	not applicable

### Other information

Bulk density:	no data available
Refraction index:	no data available
Dissociation constant:	no data available
Surface tension:	no data available
Henry's Law Constant:	no data available

## SECTION 10: Stability and reactivity

### Reactivity

This material is non-reactive under normal conditions.

### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature).

**Possibility of hazardous reactions**

No further relevant information available.

**Conditions to avoid**

No further relevant information available.

**Incompatible materials:**

No further relevant information available.

**Hazardous decomposition products**

Decomposition products in case of fire: see section 5.

**SECTION 11: Toxicological information****Information on toxicological effects****Acute effects**

*Acute oral toxicity:*

no data available

*Acute dermal toxicity:*

no data available

*Acute inhalation toxicity:*

no data available

**Irritant and corrosive effects:**

*Primary irritation to the skin:*

not applicable

*Irritation to eyes:*

not applicable

*Irritation to respiratory tract:*

not applicable

**Respiratory or skin sensitization**

In case of skin contact: not sensitizing

In case of inhalation: not sensitizing

**STOT-single exposure**

not applicable

**STOT-repeated exposure**

not applicable

**CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)****Carcinogenicity**

No indication of human carcinogenicity.

**Germ cell mutagenicity**

No indications of human germ cell mutagenicity exist.

**Reproductive toxicity**

No indications of human reproductive toxicity exist.

**Aspiration hazard**

not applicable

**Other adverse effects**

no data available

**SECTION 12: Ecological information****Ecotoxicity****Fish toxicity:**

no data available

**Daphnia toxicity:**

no data available

**Algae toxicity:**

no data available

**Bacteria toxicity:**

no data available

**Persistence and degradability**

no data available

**Bioaccumulative potential**

Partition coefficient: n-octanol/water: no data available

**Mobility in soil:**

no data available

**Other adverse effects**

no data available

**SECTION 13: Disposal considerations****Waste treatment methods****Appropriate disposal / Product**

Dispose according to legislation. Consult the appropriate local waste disposal expert about waste disposal.

**Appropriate disposal / Package**

Dispose according to legislation. Handle contaminated packages in the same way as the substance itself.

**SECTION 14: Transport information****Land transport (DOT)**

No dangerous good in sense of this transport regulation.

**Sea transport (IMDG)**

No dangerous good in sense of this transport regulation.

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code  
not relevant

**Air transport (ICAO-TI / IATA-DGR)**

No dangerous good in sense of this transport regulation.

<b>SECTION 15: Regulatory information</b>
---

**Safety, health and environmental regulations/legislation specific for the substance or mixture**

**National regulations**

**Toxic Substances Control Act (TSCA)**

Listed

**OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)**

Not listed.

**SARA 313 Components**

Not listed.

**US State Regulations**

**Massachusetts Right To Know Components**

Listed

**Pennsylvania Right To Know Components**

Listed

**New Jersey Right To Know Components**

Listed

**California Prop. 65 Components**

Not listed.



## SECTION 16: Other information

### Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists  
DOT - Department of Transportation  
IARC - International Agency for Research on Cancer  
IATA-DGR - International Air Transport Association-Dangerous Goods Regulations  
ICAO-TI - International Civil Aviation Organization-Technical Instructions  
IMDG - International Maritime Code for Dangerous Goods  
LTV - Long Term Value  
NIOSH - National Institute for Occupational Safety and Health  
NTP - National Toxicology Program  
OSHA - Occupational Safety & Health Administration  
PBT - Persistent, Bioaccumulative and Toxic  
PEL - Permissible Exposure Limit  
STV - Short Term Value  
SVHC - Substances of Very High Concern  
TDG - Transport of Dangerous Goods  
TLV - Threshold Limit Value  
vPvB - very Persistent, very Bioaccumulative

### Key literature references and sources for data

This Safety Data Sheet has been prepared based on information available for public as TOXNET information, European Chemicals Agency (ECHA) substance dossier, papers from international cancer research institutes (IARC Monographs), U.S. National Toxicology Program data, U.S. Agency for Toxic Substances and Disease Control (ATSDR), PubChem websites and SDS from our raw material manufacturers.

Revision date	Version	Print date
06.09.2023	1.1	2023-09-06

### Additional information

Indication of changes: Section 2

If you need an explanation of the change, contact the supplier (SDS@avantorsciences.com).

### Disclaimer

The information provided in this Safety Data Sheet (SDS) was prepared based on data believed to be accurate as of the date of this SDS. TO THE GREATEST EXTENT PERMITTED BY LAW, AVANTOR PERFORMANCE MATERIALS ("AVANTOR") EXPRESSLY DISCLAIMS ANY AND ALL REPRESENTATIONS AND WARRANTIES REGARDING THE INFORMATION CONTAINED HEREIN INCLUDING, WITHOUT LIMITATION, AS TO ACCURACY, COMPLETENESS, FITNESS FOR PURPOSE OR USE, MERCHANTABILITY, NON-INFRINGEMENT, PERFORMANCE, SAFETY, SUITABILITY AND STABILITY. This SDS is intended as a guide to the appropriate use, handling, storage and disposal of the product to which it relates by properly trained personnel, and is not intended to be comprehensive. Users of Avantor's products are advised to perform their own tests and to exercise their own judgment to determine the safety, suitability and appropriate use, handling, storage and disposal of each product and product combination for their own purposes and uses. TO THE GREATEST EXTENT PERMITTED BY LAW, AVANTOR DISCLAIMS LIABILITY FOR, AND BY USING AVANTOR'S PRODUCTS PURCHASER AGREES THAT UNDER NO CIRCUMSTANCES SHALL AVANTOR BE LIABLE FOR, SPECIAL, INDIRECT, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES OF ANY TYPE OR KIND, INCLUDING WITHOUT LIMITATION, FOR LOSS OF PROFITS, REPUTATIONAL DAMAGE, PRODUCT RECALL OR BUSINESS INTERRUPTION.

